## Microscopi Service Rules

### Remember:

- All equipment at the SCC-PCB facility are in self-service mode.
- Before using any equipment, you must read the manual available in the room.
- Not following these recommendations will result in a misuse of the equipment / room

### 1. The service:

The microscopy room is located at EMS1BC36. This service is intended for the labs installed at the PCB with a the correspondent SCC-PCB fee. Room's access is restricted by card. To use the equipment, prior training is required by the person in charge of the service. To ensure proper operation SCC-PCB staff check the facilities daily.

### 2. Use of Microscope:

#### 2.1. New users

You must register in the service by filling out the form on the PCB intranet: *Microscopy user registration Form*. Once your request is received, SCC-PCB Staff will contact you to schedule a training session. The training takes place from Monday to Friday at 3 p.m.

Once the training is done, your card will get access granted for the room.

#### 2.2. Microscopy Equipment Reservation:

- Is mandatory to book the equipment by using the SCC-PCB online booking system.
- On working days from 9am to 6pm, bookings of more than 3 hours are not allowed.
- Booking that will not be used must be cancelled. If you can not delete it, please contact by e-mail to [scc-pcb@pcb.ub.es](mailto:scc-pcb@pcb.ub.es).
- Remember that your username and password are personal and non-transferable.

### 3. Use of Microscopy

- Never touch the microscopes with gloves.
- Never force any part of the microscope: revolver, plate, objectives, eyepiece, condenser, etc.
- At the end of the session, the eyepieces and lenses must be cleaned with KIMTECH™ wipes soaked with a little 70% ethanol.
- It is mandatory to register the use of the microscope in the equipment logbook indicating: use, fluorescence use, hours of use and filters used.
- If the epifluorescence light source is cold, protect the microscope with its cover.

### 4. Safety rules

#### **RISK OF EXPLOSION**

Mercury lamps can explode releasing highly toxic mercury vapor into the environment. This risk occurs when the lamp has no time to cool. To avoid these accident, follow the instructions below when handling epifluorescence light source:

- Before turning on the fluorescence light source**, check how long it has been used before you do. Please, wait 30 minutes after the last use. You can check the time on the logbook or on the online booking system.
- When you are finished using the epifluorescence light source**, check on the logbook or on the online booking system if any user has booked it for the next 30 minutes. If not, you can close the epifluorescence light source

### 5. Data storage:

To prevent computer viruses or computer problems in the microscope PC, the USB ports are disabled. **USB can not be used to transfer data.**

PCB user must transfer data through the PCB's CABALLA Server:

- Images will be saved on the server for a maximum of 1 week. After this time or if the server is full, they are automatically deleted.
- To access the server:

<p><b>Windows user:</b></p> <p>In the Windows explorer search: <a href="https://caballa.adm.pcb.ub.es">\\caballa.adm.pcb.ub.es</a></p> <p>Click on the “microscopy” folder</p> <p>Enter username: <b>adm\username</b>, and enter the password</p> <p><b>You can never select the option: Remember password</b></p>	<p><b>Mac user:</b></p> <p>Finder/Go/Connect on the Server, type in the address bar: <b>smb: \\caballa.adm.pcb.ub.es</b></p> <p>Enter username <b>adm\username</b>, and enter the password</p> <p>Click on “microscopy” folder</p> <p><b>You can never select the option: Remember password</b></p>
--	---

We are at your disposal for any questions or clarifications.

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